


# OptiPrep™ Application Sheet C11

## Isolation of human polymorphonuclear leukocytes (granulocytes) from a leukocyte rich plasma in a discontinuous iodixanol gradient

- ◆ OptiPrep™ is a 60% (w/v) solution of iodixanol in water, density = 1.32 g/ml
- ◆ **Axis-Shield Mini-Review (MC01)** “Purification of mononuclear cells, monocytes and polymorphonuclear leukocytes – a methodological review” compares all of the currently available methodologies
- ◆ **Axis-Shield Mini-Review (MC02)** “Purification of mononuclear cells, monocytes and polymorphonuclear leukocytes – a bibliographical review” provides a comprehensive bibliography of all the published papers reporting the use of OptiPrep™
- ◆ To access **MC01 and MC02** return to the initial list of Folders and select “**Mini-Reviews**”.
- ◆ To access other Application Sheets referred to in the text return to the Cell Index; key Ctrl “F” and type the C-Number in the Find Box



**Figure 1:** Density of human blood cells

erythrocytes; creating a continuous gradient in which the MCs and PMNs band according to their buoyant densities. The theory and methodology of this strategy is fully described in the Polymorphprep™ Application Sheet.

There are two alternative strategies that involve the use of a simple density barrier ( $\rho = 1.077$  g/ml), either Lymphoprep™ or a Nycodenz® solution. (1) Whole blood is layered over the density barrier; the MCs band at the interface while the erythrocytes and PMNs sediment through the barrier to form a pellet. The latter is then recovered and the erythrocytes selectively lysed in isotonic  $\text{NH}_4\text{Cl}$  or ice-cold water to leave a pure PMN fraction. (2) A leukocyte-rich plasma (LRP), prepared from whole blood by polysucrose-aggregation of the erythrocytes, is layered over the density barrier, which then separates the MCs and PMNs.

- ◆ To access Polymorphprep™ or Lymphoprep™ Application Sheets go to the Axis-Shield Home Page and click on “Products”.

Pelleting of PMNs however can easily cause aggregation and activation of neutrophils and it is strongly recommended that if an LRP is used, a dense cushion (1.090-1.095 g/ml) be placed beneath the 1.077 g/ml layer in order to band the PMNs rather than pellet them. The two density gradient solutions can be easily prepared from OptiPrep™ and it is this method that is described in the Application Sheet.

### 1. Background

With the exception of basophils, the polymorphonuclear leukocytes (PMNs) or granulocytes from human peripheral blood have densities predominantly above 1.080 g/ml, while mononuclear cells (MCs) have densities below 1.077 g/ml (see Figure 1). Since the density of erythrocytes significantly overlaps that of the denser neutrophils there is only one means by which PMNs may be isolated from whole human blood using a single step method and that is to use the Axis-Shield medium called Polymorphprep™. It contains diatrizoate and a polysaccharide [1] and has a high density (1.113 g/ml) and osmolality (445 mOsm). The high osmolality causes loss of water from the

- ◆ Preparation of the leukocyte-rich plasma is achieved by adjusting the blood to 0.6% (w/v) polysucrose and allowing the aggregated erythrocytes to sediment at 1 g. There are however several popular alternatives that are described in Notes 1 and 2.

## 2. Solutions required

- OptiPrep™ (shake gently before use).
- Diluent: 0.85% (w/v) NaCl, 1 mM EDTA, 20 mM HEPES-NaOH, pH 7.4
- Polysucrose: 6% (w/v) polysucrose ( $M_r = 400-500 \times 10^3$ ) in 0.85% (w/v) NaCl
- Lysis buffer: 0.83% (w/v)  $\text{NH}_4\text{Cl}$ , 10 mM HEPES-NaOH, pH 7.0
- 1.8 (w/v) NaCl, 20 mM Hepes-NaOH, pH 7.4

- ◆ Strategies for preparing density solutions for mammalian cells are described in **Application Sheet C01**.

## 3. Protocol

- To 9 vol. of freshly drawn blood (containing 2 mM EDTA as anticoagulant) add 1 vol. of Solution C (see Notes 1 and 2).
- Allow the aggregated erythrocytes to settle to the bottom (20-40 min at room temperature); then aspirate the entire supernatant.
- Prepare the following density solutions from OptiPrep™ and Solution B (respectively): 1.077 g/ml, 5 vol. + 17 vol. and **EITHER** 1.090 g/ml, 8 vol. + 22 vol.; **OR** 1.095 g/ml, 17 vol. + 43 vol. (see Notes 3 and 4).
- Underlayer 5 ml of LRP with 2.5-3.0 ml of 1.077 g/ml solution and the same volume of **EITHER** 1.090 g/ml **OR** 1.095 g/ml (see Fig. 2 and Note 5).
- Centrifuge at 18-22°C for 25 min at 800 g.
- Harvest the PMNs from the lower interface and the mononuclear cells from the upper interface (see Fig. 2).
- Dilute the PMN suspension with an equal volume of Solution B and collect the PMNs by centrifugation at 250-350 g for 10 min.
- Resuspend the pellet in a suitable medium for analysis.
- To remove residual erythrocyte contamination of the PMNs, resuspend the cell pellet in 3 ml of Solution D and incubate at 37°C for 7 min **OR** resuspend the PMNs in 3 ml ice-cold distilled water, then after 30 sec add an equal volume of Solution E.
- Harvest the PMNs by centrifugation and resuspend in a suitable medium (see Note 6).


Keep the following stock solutions at 4°C  
 100 mM HEPES (free acid) 2.38 g per 100 ml water  
 100 mM EDTA( $\text{Na}_2 \cdot 2\text{H}_2\text{O}$ ) 3.72 g per 100 ml water

Solution B: Dissolve 0.85 g NaCl in 50 ml water, add 20 ml and 1 ml of HEPES and EDTA stock solutions respectively; adjust to pH 7.4 with 1 M NaOH; make up to 100 ml.

Solution C: Add the polysucrose **slowly** to the rapidly stirred NaCl solution.

Solution D: Dissolve 0.83 g  $\text{NH}_4\text{Cl}$  in 50 ml water; add 10 ml of HEPES stock solution, adjust to pH 7.0 with 1 M NaOH and make up to 100ml.

Solution E: Dissolve 1.8 g NaCl in 50 ml water, add 20 ml of HEPES stock solution; adjust to pH 7.4 with 1 M NaOH; make up to 100 ml



**Figure 2** Separation of polymorphs from mononuclear cells from a leukocyte-rich plasma.

#### 4. Notes

1. If exposure of the cells to polysucrose is deemed undesirable (there is evidence that this macromolecule can adsorb to the surface of leukocytes) then centrifuge the blood at  $200 g_{av}$  for 15-20 min at 18-22°C and harvest the buffy coat in the plasma supernatant from the top of the packed erythrocytes and use this in Step 4.
2. Methylcellulose can be added to the blood as an alternative to polysucrose for erythrocyte aggregation. A variant of the method for removing the erythrocytes, adapted from Boyum [2], is to layer the blood over 12% (w/v) iodixanol in 130 mM NaCl, containing 1.66% methylcellulose. After standing at 1g the aggregated erythrocytes sediment to the bottom of the tube [3,4].
3. If the density of this cushion is 1.090 g/ml a small percentage of the neutrophils and most of the eosinophils will sediment through this layer. If a density of 1.095 g/ml is chosen, virtually all of the PMNs will be retained by the high-density barrier. On the other hand, fewer of the residual erythrocytes in the LRP will contaminate the PMN band using the lower density cushion. See section 5 for some examples.
4. Occasionally the low-density 1.077 g/ml iodixanol layer is replaced with a routine commercial peripheral blood mononuclear (PBMC) isolation medium such as Lymphoprep™. There may be some merit in the use of this solution as one of its components, 5.7% (w/v) polysucrose, may cause the residual erythrocytes to aggregate further and thus assist their sedimentation into a pellet. On the other hand this concentration of polysucrose is almost ten times higher than that to which the blood is adjusted in Step 2 of the protocol. See Section 5 for some examples.
5. For more information about preparing discontinuous gradients see [Application Sheet C02](#).
6. In a recent study of twelve healthy males [5], the two-layer iodixanol gradient consistently gave a yield of  $\geq 95\%$ .

#### 5. Methodological review

Some of the gradient and centrifugation conditions described in published papers reporting the use of the method described in this Application Sheet to purify PMNs from primate peripheral blood are summarized in Table 1.

**Table 1** Primate PMN isolation using OptiPrep™

Species	Gradient and centrifugation comments	Ref. #
Human	Iodixanol solutions of 1.077 + 1.090 g/ml, 1200 g for 25 min	6
	Iodixanol solutions of 1.077 + 1.090 g/ml, 800 g for 30 min	7
	Iodixanol solutions of 1.077 + 1.090 g/ml, 600 g for 20 min	8
	PBMC isolation medium + iodixanol solution of 1.095 g/ml, 400/500 g for 25/20 min	9-16
Monkey	Iodixanol solutions of 1.077 + 1.095 g/ml	17

The methodology has been used, with minor variations, for the isolation of PMNs from bovine, guinea pig, mouse, rabbit and rat sources. For references see Mini-Review MC01.

#### 6. References

1. Ferrante, A. and Thong, Y.H. (1980) *Optimal conditions for simultaneous purification of mononuclear and polymorphonuclear leucocytes from human blood by the Hypaque-Ficoll method* J. Immunol. Meth., **36**, 109-117
2. Boyum, A. (1968) *Isolation of mononuclear cells and granulocytes from human blood: Isolation of mononuclear cells by one centrifugation and of granulocytes by combining centrifugation and sedimentation at 1g* Scand. J. Clin. Lab. Invest., **21**, 1-89
3. Niggli, V. (2003) *Microtubule-disruption-induced and chemotactic-peptide-induced migration of human neutrophils: implications for differential sets of signaling pathways* J. Cell Sci., **116**, 813-822
4. Dehghani Zadeh, A., Seveau, S., Halbwachs-Mecarelli, L. and Keller, H.U. (2003) *Chemotactically-induced redistribution of CD43 as related to polarity and locomotion of human polymorphonuclear leucocytes* Biol. Cell, **95**, 265-273
5. Radom-Aizik, S., Zalvidar, Jr., F., Leu, S-Y., Galasetti, P. and Cooper, D.M. (2008) *Effects of 30 min of aerobic exercise on gene expression in human neutrophils* J. Appl. Physiol., **104**, 236-243
6. Yuan, Z-N, Tolo, K., Schenck, K. and Helgeland, K. (1999) *Increased levels of soluble Fcγ receptor III in gingival fluid periodontal lesions* Oral Microbiol. Immunol., **14**, 172-175

7. Hudgens, J., Langkamp-Henken, B., Stechmiller, J.K., Herrlinger-Garica, K.A. and Nieves, C. (2005) *Immune function is impaired with a mini nutritional assessment score indicative of malnutrition in nursing home elders with pressure ulcers* J. Parenter. Enteral Nutr. **28**, 416-422
8. Chiu, H-C., Liang, J-S., Wang, J-S. and Lu, J-F. (2006) *Mutational analyses of Taiwanese kindred with X-linked adrenoleukodystrophy* Pediatr. Neurol., **35**, 250-256
9. Channon, J.Y., Seguin, R.M. and Kasper, L. (2000) *Differential infectivity and division of Toxoplasma gondii in human peripheral blood leukocytes* Infect. Immun., **68**, 4822-4826
10. Channon, J.Y., Miselis, K.A., Minns, L.A, Dutta, C. and Kasper, L.H. (2002) *Toxoplasma gondii induces granulocyte colony-stimulating factor and granulocyte-macrophage colony-stimulating factor secretion by human fibroblasts: implications for neutrophil apoptosis* Infect. Immun., **70**, 6048-6057
11. Shen, L., Fahey, J.V., Hussey, S.B., Asin, S.N., Wira, C.R. and Fanger, M.W. (2004) *Synergy between IL-8 and GM-CSF in reproductive tract epithelial cell secretions promotes enhanced neutrophil chemotaxis* Cell. Immunol., **230**, 23-32
12. Shen, L., Smith, J.M., Shen, Z., Hussey, S.B., Wira, C.R. and Fanger, M.W. (2006) *Differential regulation of neutrophil chemotaxis to IL-8 and fMLP by GM-CSF: lack of direct effect of oestradiol* Immunology, **117**, 205-212
13. Smith, J.M., Wira, C.R., Fanger, M.W. and Shen, L. (2006) *Human fallopian tube neutrophils – a distinct phenotype from blood neutrophils* Am. J. Reprod. Immunol., **56**, 218-229
14. Shen, L., Smith, J.M., Shen, Z., Eriksson, M., Sentman, C. and Wira, C.R. (2007) *Inhibition of human neutrophil degranulation by transforming growth factor- $\beta$ 1* Clin. Exp. Immunol., **149**, 155-161
15. Smith, J.M., Shen, Z., Wira, C.R., Fanger, M.W. and Shen, L. (2007) *Effects of menstrual cycle status and gender on human neutrophil phenotype* Am. J. Reprod. Immunol., **58**, 111-119
16. Pioli, P.A., Jensen, A.L., Weaver, L.K., Amiel, E., Shen, Z., Shen, L., Wira, C.R. and Guyre, P.M. (2007) *Estradiol attenuates lipopolysaccharide-induced CXC ligand 8 production by human peripheral blood monocytes* J. Immunol., **179**, 6284-6290
17. Lau, M., Vayntrub, T., Grumet, F.C., Lowsky, R., Strober, S., Hoppe, R., Larson, M., Holm, B., Reitz, B. and Borie, D (2004) *Short tandem repeat analysis to monitor chimerism in Macaca Fascicularis* Am. J. Transplant., **4**, 1543-1548

Application Sheet C11; 7<sup>th</sup> edition, June 2016

