

OptiPrep™ Application Sheet C23

Isolation of a progenitor cell-enriched fraction methods

- ◆ OptiPrep™ is a 60% (w/v) solution of iodixanol in water, density = 1.32 g/ml
- ◆ To access other Application Sheets referred to in the text return to the Cell Index; key Ctrl “F” and type the C-Number in the Find Box
- ◆ Section 1 describes the more traditional sedimentation from a cell suspension on to a density barrier; Section 2 describes the use of the sample itself as the density barrier; Section 3 briefly describes a strategy for neural tissue progenitors

1. Sedimentation on to a density barrier

1a Introduction

The study of the proliferation and maturation of stem cells from bone marrow is a major area of cell and molecular biology research. Usually the resolving power of simple density barriers is insufficient to allow progenitor cells to be isolated in a sufficiently pure form for further analysis and culture. Nevertheless, the low-density cell fraction harvested from the barrier interface can provide a progenitor cell-enriched population that allows a more effective and economical use of immuno-magnetic beads to remove lineage cells (using a cocktail of lineage-specific monoclonal antibodies).

Although metrizamide was used for density barriers prior to 1984, the lower cell toxicity of Nycodenz®, made this the medium of choice. Mayanagi et al [1] also emphasized the lack of toxicity of Nycodenz® (compared to Percoll®) and the avoidance of positive selection with antibody-beads was important in recovering viable progenitor cells. Bertoncello et al [2,3] reported the use of a 1.085 g/ml Nycodenz® barrier to isolate low-density cells from a suspension of nucleated bone marrow cells; 3 ml containing 10^8 cells were layered over 5 ml of the density barrier and centrifuged at 100 g for 30 min, at 4°C and the progenitor cell-enriched fraction harvested from the interface. A barrier of the same density was used by Erlich et al [4] and Sitnicka et al [5]. Lower density barriers, 1.080 g/ml [6], 1.077 g/ml [7-10] and 1.068 g/ml [11] have also been used and frequently the centrifugation is carried out at room temperature rather than at 4°C. In some cases, a shorter centrifugation time (20 min) has been employed at 1000 g [12]; in others a g-force of 400 g was used [5,6]. The 1.077 g/ml density barrier has been widely used since ready-made solutions of this density have been available as NycoPrep™ 1.077 and NycoPrep™ 1.077A, the latter being the more popular for murine samples (see also Section 3). The density barrier approach has been adapted to OptiPrep™; several publications report the use of a $\rho = 1.08$ g/ml iodixanol barrier [13-17]. As with the Nycodenz® methodology however, there are a number of variations for the density of the iodixanol barrier: 1.077 g/ml [18-21] and 1.074 g/ml [22].

Progenitor-enriched fractions have also been isolated from human blood [23] and cord blood [24] using NycoPrep™ 1.077; those from murine blood [25,26] and fetal liver [27-30] on NycoPrep™ 1.077A. Again the g-force is rather variable – 400-1000 g. Using OptiPrep™, progenitor cells have been enriched from blood (1.072 g/ml) [31] and thymus (1.05 g/ml) [32-34].

- ◆ Important note: As many Nycodenz® gradient solutions were routinely used as NycoPrep™ 1.077A or diluted from NycoPrep™ 1.15, neither of which are now commercially available, the following methodology is based on the use of OptiPrep™ (a sterile solution of 60% w/v iodixanol) to avoid making up solutions from Nycodenz® powder.
- ◆ Modified density barriers, in which the density of the sample is raised before layering on the barrier, or the use of discontinuous gradients may provide higher enrichments of progenitor cells than a simple density barrier. For more information on the advantages of using such techniques see “Hepatic cells”, Application Sheet C47 in index.

1b. Solutions required (see Note 1)

- A. OptiPrep™ (shake the bottle gently before use)
- B. Iscove's modified Dulbecco's medium Hank's (IMDM)
- C. Fetal calf serum (FCS)
- D. Phosphate buffered saline (PBS) **OR**
- E. Tricine buffered saline

Keep Tricine as 100 mM stock solution at 4°C; 1.79g per 100 ml water.

Solution E: Dissolve 0.85g NaCl in 50 ml water; add 10 ml of Tricine stock solution; adjust to pH 7.0 with 1 M NaOH and make up to 100 ml.

1c. Protocol

Carry out all operations at room temperature (see Note 2)

1. Harvest bone marrow material in Solution B containing 2% FCS.

2. Dissociate the bone marrow cells by passing three times through a 25-gauge needle.

Table 1 Density of iodixanol solutions prepared by dilution of OptiPrep with Solutions D or E

Density (g/ml)	OptiPrep	Solution D/E
1.050	8.5 vol.	51.5 vol.
1.063	11.0 vol.	49.0 vol.
1.068	12.0 vol.	48.0 vol.
1.077	13.5 vol.	46.5 vol.
1.080	14.0 vol.	46.0 vol.
1.085	15.0 vol.	45.0 vol.
1.090	16.0 vol.	44.0 vol.

3. Filter the suspension through 200 µm nylon mesh.

4. To prepare a density barrier solution of approx 395 mOsm, dilute OptiPrep™ directly with Solution D or E (with or without 5% FCS) to give a solution of the chosen density (see Table 1 and Notes 3 and 4).

5. To prepare a density barrier solution of 1.077 g/ml and osmolality approx 265 mOsm (equivalent to Nycoprep 1.077A) start by diluting Solution E with water (volume ratio 2.5:0.5 respectively); this solution has an osmolality of approx 242 mOsm. Dilute OptiPrep™ with the 242 mOsm solution using a volume ratio of 2.7:9.3 respectively.

6. Layer 5 ml of the bone marrow cell suspension (10^7 cells/ml) over 3 ml of the chosen density barrier (see Note 5).

7. Centrifuge at 450 g for 20 min; allow the rotor to decelerate without the brake (see Notes 6 and 7).

8. Collect the cells from the interface; dilute with an equal volume of Solution D or E (+ 5% FCS); wash twice and resuspend in the same medium prior to lineage committed cell depletion.

1d. Notes

1. To dilute the OptiPrep™ use any suitable medium, two examples are given as Solutions D or E; culture medium may also be used. As long as the solution is isoosmotic, then the dilutions will also be isoosmotic.
2. Some papers report the use of centrifugation temperatures of 4°C.
3. Step 5 describes the preparation of an alternative reduced osmolality OptiPrep™ density solution.
4. Dilutions of OptiPrep™ with any isotonic medium (including culture medium) will give the same densities as those shown in Table 1. These diluents all have a density of approx 1.006 g/ml. Inclusion of FCS will however modulate the density slightly upwards. FCS has a density of approx 1.032 g/ml, so any balanced salt solution containing 10% FCS for example will have a density of 1.009 g/ml. The density of the OptiPrep™/saline mixtures given in Table 1 will be increased by 0.002-0.003 g/ml. If 5% FCS is used the density difference is proportionately less. For more information on preparing gradient solutions see [Application Sheet C01](#).
5. Since the megakaryocytic progenitor cells are less dense than the majority of other cell types, it may be worth considering using a flotation strategy in which the bone marrow cell suspension is adjusted to the barrier density (e.g. 1.085 g/ml) and the progenitor cells harvested from the top of the solution. This approach has been used for fetal liver cells [35].

6. The published g-force, time and temperature values are rather variable. What effect lower or higher values might have on the separation can only be determined by experimentation.
7. Using the brake may create a vortex in the liquid and cause loss of definition of the interfacial band and even contamination from denser cells.

2. Progenitor cells from blood using a mixer strategy

This was first introduced for the isolation of mononuclear cells from human blood and later rat and mouse blood (see **“Mononuclear cells” Application Sheets C04, C06 and C07 respectively in index**). It has also been applied to the isolation of progenitor cells. The crude cell suspension is PBS + 5% FCS (30 ml) was mixed with 10.1 ml of OptiPrep™ and 4 ml of water. This adjusts the cell suspension medium to approx 1.077 g/ml. A small volume of buffered saline is placed on top to prevent the low-density cells from banding at an air-liquid interface during centrifugation [36-41]. The strategy avoids the build up and possible aggregation of cells at the interface between the sample and the density barrier in the standard barrier sedimentation strategy (see Section 1). It allows the denser cells to sediment while the low-density cells float. The 4 ml of water used in the preparation of the sample could certainly be replaced with saline if required, without affecting the separation.

- ◆ In a symposium abstract Hu et al [42] also reported a new OptiPrep™-based method for the purification of blood-derived mesenchymal stem cells.

3. Discontinuous gradients

A novel method [43,44] used gradient solutions produced by mixing OptiPrep™ with NycoPrep 1.077A. The latter is now commercially unavailable but a solution of the same density and reduced osmolality can be produced from OptiPrep™ by dilution with a reduced osmolality saline solution see **“Mononuclear cells” Application Sheet C43 in index for more information**) Solutions can be prepared as follows: use 0.7% (w/v) NaCl, 10 mM Tricine, pH 7.0 to dilute OptiPrep™ to produce solutions of 1.050, 1.080, 1.090 g/ml (approx. iodixanol concentrations if 8.4, 14 and 16% w/v). The gradient is centrifuged at 400 g for 15 min to isolate a low density fraction.

The iodixanol density gradient devised by Brewer et al [45] for the isolation of hippocampal motoneurons is also used as a gradient for purifying neural progenitor cells [46-50]. **The method is described in detail in “Neural Cells” Application Sheet C29 in index.** Abbosh et al [50] reported a two-step gradient of 1.058 and 1.11 g/ml for postnatal hippocampal precursors; He and Shen [52] and Nunan et al [53] also used a discontinuous iodixanol gradient for the isolation of glial progenitor cells, but gave no details.

The oval cells of liver, which are considered to be the equivalent of progenitor cells in the liver, may potentially be propagated *in vitro* and used restoratively in some liver diseases. Three layer gradients of 13%, 16 and 18% (w/v) iodixanol (OptiPrep™ diluted with HBSS, supplemented with 0.2% BSA (10 ml, 10 ml and 5 ml respectively) are overlaid with the non-parenchymal cells (see Application Sheet C24) in 11% iodixanol and centrifuged at 6,500 g for 30 min [54]. The oval cells are enriched at the 11/13% interface.

4. References

1. Mayanagi, T., Kurosawa, R., Ohnuma, K., Ueyama, A., Ito, K. and Takahashi, J. (2003) *Purification of mouse primordial germ cells by Nycodenz* Reproduction, **125**, 667-675
2. Bertocello, I., Bartelmez, S.H., Bradley, T.R. and Hodgson, G.S. (1987) *Increased Qa-m7 antigen expression is characteristics of primitive hemopoietic progenitors in regenerating marrow* J. Immunol., **139**, 1096-1103
3. Bertocello, I., Bradley, T.R., Hodgson, G.S. and Dunlop, J.M. (1991) *The resolution, enrichment, and organization of normal bone marrow high proliferative potential colony-forming cell subsets on the basis of rhodamine-123 fluorescence* Exp. Hematol., **19**, 174-178
4. Erlich, S., Miranda, S.R.P., Visser, J.W.M., Dagan, A., Gatt, S. and Schuchman, E.H. (1999) *Fluorescence-based selection of gene-corrected hematopoietic stem and progenitor cells from acid sphingomyelinase-deficient mice: implications for Niemann-Pick disease gene therapy and the development of improved stem cell gene transfer procedures* Blood, **93**, 80-86

5. Sitnicka, E., Ruscetti, F.W., Priestley, G.V., Wolf, N.S. and Bartelmez, S.H. (1996) *The effect of thrombopoietin on the proliferation and differentiation of murine hematopoietic stem cells* Blood, **88**, 82-88
6. Hodohara, K., Fujii, N., Yamamoto, N. and Kaushansky, K. (2000) *Stromal cell-derived factor-1 (SD-1) acts together with thrombopoietin to enhance the development of megakaryocytic progenitor cells (CFU-MK)* Blood, **95**, 769-775
7. Morgan, B., Sun, L., Avitahl, N., Andrikopoulos, K., Ikeda, T., Gonzales, E., Wu, P., Neben, S. and Georgopoulos, K. (1997) *Aiolos, a lymphoid restricted transcription factor that interacts with Ikaros to regulate lymphocyte differentiation* EMBO J., **16**, 2004-2013
8. Nibbs, R.J.B., Wylie, S.M., Pragnell, I.B. and Graham, G.J. (1997) *Cloning and characterization of a novel murine β chemokine receptor, D6* J. Biol. Chem., **272**, 12495-12504
9. Tsai, S., Fero, J. and Bartelmez, S. (2000) *Mouse jagged2 is differentially expressed in hematopoietic progenitors and endothelial cells and promotes the survival and proliferation of hematopoietic progenitors by direct cell-to-cell contact* Blood, **96**, 95-957
10. Weich, N.S., Fitzgerald, M., Wang, A., Calveti, J., Yetz-Aldape, J., Neben, S. and Turner, K.J. (2000) *Recombinant human interleukin-11 synergizes with steel factor and interleukin-3 to promote directly the early stages of murine megakaryocyte development in vitro* Blood, **95**, 503-509
11. Berkovic, D., Bensch, M., Bertram, J., Wille, T., Haase, D., Binder, C. and Fleer, E.A.M. (2001) *Effects of hexadecylphosphocholine on thrombocytopoiesis* Eur. J. Cancer, **37**, 503-511
12. Nichogiannopoulou, A., Trevisan, M., Neben, S., Friedrich, C. and Georgopoulos, K. (1999) *Defects in hemopoietic stem cell activity in Ikaros mutant mice* J. Exp. Med., **190**, 1201-1213
13. Rojnuckarin, P., Drachman, J.G. and Kaushansky, K. (1999) *Thrombopoietin-induced activation of the mitogen-activated protein kinase (MAPK) pathway in normal megakaryocytes: role in endomitosis* Blood, **94**, 1273-1282
14. Hodohara, K., Fujii, N., Yamamoto, N. and Kaushansky, K. (2000) *Stromal cell-derived factor-1 (SDF-1) acts together with thrombopoietin to enhance the development of megakaryocytic progenitor cells (CFU-MK)* Blood, **95**, 769-775
15. Geddis, A.E., Fox, N.E. and Kaushansky, K. (2001) *Phosphatidylinositol 3-kinase is necessary but not sufficient for thrombopoietin-induced proliferation in engineered Mpl-bearing cell lines as well as in primary megakaryocytic progenitors* J. Biol. Chem. **276**, 34473-34479
16. Wang, Q., Miyakawa, Y., Fox, N. and Kaushansky, K. (2000) *Interferon- α directly represses megakaryopoiesis by inhibiting thrombopoietin-induced signaling through induction of SOCS-1* Blood, **96**, 2093-2099
17. Carow, C.E., Fox, N.E. and Kaushansky, K. (2001) *Kinetics of endomitosis in primary murine megakaryocytes* J. Cell. Physiol., **188**, 291-303
18. Whetton, A.D., Lu, Y., Pierce, A., Carney, L. and Spooncer, E. (2003) *Lysophospholipids synergistically promote primitive hematopoietic cell chemotaxis via a mechanism involving Vav1* Blood, **102**, 2798-2802
19. Evans, C.A., Tonge, R., Blinco, D., Pierce, A., Shaw, J., Lu, Y., Hanzah, H.G., Gray, A., Downes, C.P., Gaskell, S.J., Spooncer, E. and Whetton, A.D. (2004) *Comparative proteomics of primitive hematopoietic cell populations reveals differences in expression of proteins regulating motility* Blood, **103**, 3751-3759
20. Unwin, R.D., Smith, D.L., Blinco, D., Wilson, C.L., Miller, C.J., Evans, C.A., Jaworska, E., Baldwin, S.A., Barnes, K., Pierce, A., Spooncer, E. and Whetton, A.D. (2006) *Quantitative proteomics reveals posttranslational control as a regulatory factor in primary hematopoietic stem cells* Blood, **107**, 4687-4694
21. Dooner, M.S., Aliotta, J.M., Pimental, J., Dooner, G.J., Abedi, M., Colvin, G., Liu, Q., Weier, H-U., Johnson, K. and Quesenberry, P.J. (2008) *Conversion potential of marrow cells into lung cells fluctuates with cytokine-induced cell cycle* Stem Cells Devel., **17**, 207-219
22. Broudy, V.C. and Lin, N.L. (2004) *AMG531 stimulates megakaryopoiesis in vitro by binding to Mpl* Cytokine, **25**, 52-60
23. Hepburn, M.D., Nagesh, K., Heppleston, A.D., Cachia, P.G. and Pippard, M.J. (2001) *Timing of the appearance of multipotential and committed haemopoietic progenitors in peripheral blood after mobilization in patients with lymphoma* Clin. Lab. Hematol., **23**, 119-124
24. Zeng, F., Chen, M-J., Huang, W-Y., Yan, J-B., Xiao, Y-P., Gong, Z-J., Ren, Z-R. and Huang, S-Z. (2005) *In utero transplantation of human hematopoietic stem cells into fetal goats under B-type ultrasonographic scan: an experimental model for the study of potential prenatal therapy* Eur. J. Obstet. Gynecol. Reprod. Biol., **118**, 170-173
25. Meister, B., Maurer, H., Simma, B., Kern, H., Ulmer, H., Hittmair, A. and Fink, F-M. (1997) *The effect of recombinant human erythropoietin on circulating hematopoietic progenitor cells in anemic premature infants* Stem Cells, **15**, 359-363
26. Yamamoto, Y., Yasumizu, R., Amou, Y., Watanabe, N., Nishio, N., Toki, J., Fukuhara, S. and Ikehara, S. (1996) *Characterization of peripheral blood stem cells in mice* Blood, **88**, 445-454
27. Kurata, H., Mancini, G.C., Alespeiti, G., Migliaccio, A.R. and Migliaccio, G. (1998) *Stem cell factor induces proliferation and differentiation of fetal progenitor cells in the mouse* Br. J. Hematol., **101**, 676-687
28. Barcena, A., Muench, M.O., Song, K.S., Ohkuba, T. and Harrison, M.R. (1999) *Role of CD95/Fas and its ligand in the regulation of the growth of human CD34⁺CD38⁻ fetal liver cells* Exp. Hematol., **27**, 1428-1439
29. Huie, M.A., Cheung, M-C., Muench, M.O., Becerril, B., Kan, Y.W. and Marks, J.D. (2001) *Antibodies to human fetal erythroid cells from a non-immune phage antibody library* Proc. Natl. Acad. Sci. USA, **98**, 2682-2687
30. Shao, J., Stapleton, P.L., Lin, Y.S. and Gallagher, E.P. (2007) *Cytochrome P450 and glutathione S-transferase mRNA expression in human fetal liver hematopoietic stem cells* Drug Metab. Dispos., **35**, 168-175
31. Porat, Y., Porozov, S., Belkin, D., Shimoni, D., Fisher, Y., Belleli, A., Czeiger, D., Silverman, W.F., Belkin, M., Battler, A., Fulga, V. and Savion, N. (2006) *Isolation of an adult blood-derived progenitor cell population capable of differentiation into angiogenic, myocardial and neural lineages* Br. J. Haematol., **135**, 703-714

32. Porritt, H.E., Rumfelt, L.L., Tabrizifard, S., Schmitt, T.M., Zuniga-Pflucker, J.C. and Petrie, H.T. (2004) *Heterogeneity among DN1 prothymocytes reveals multiple progenitors with different capacities to generate T cell and non-T cell lineages* Immunity, **20**, 735-745
33. Tabrizifard, S., Oлару, A., Plotkin, J., Fallahi-Sichani, M., Livak, F. and Petrie, H.T. (2004) *Analysis of transcription factor expression during discrete stages of postnatal thymocyte differentiation* J. Immunol., **173**, 1094-1102
34. Umland, O., Mwangi, W.N., Anderson, B.M., Walker, J.C. and Petrie, H.T. (2007) *The blood contains multiple distinct progenitor populations with clonogenic B and T lineage potential* J. Immunol., **178**, 4147-4152
35. Muench, M., Cupp, J., Polakoff, J. and Roncarolo, M.G. (1994) *Expression of CD33, CD38, and HLA-DR on CD34⁺ human fetal liver progenitors with a high proliferative potential* Blood, **83**, 3170-3181
36. Moore, B.E., Colvin, G.A., Dooner, M.S. and Quesenberry, P.J. (2005) *Lineage-negative bone marrow cells travel bidirectionally in the olfactory migratory stream but maintain hematopoietic phenotype* J. Cell. Physiol., **202**, 147-152
37. Colvin, G.A., Dooner, M.S., Dooner, G.J., Sanchez-Guijo, F.M., Demers, D.A., Abedi, M., Ramanathan, M., Chung, S., Pascual, S. and Quesenberry, P.J. (2007) *Stem cell continuum: Directed differentiation hotspots* Expt. Hematol., **35**, 96-107
38. Dooner, M.S., Aliotta, J.M., Pimental, J., Dooner, G.J., Abedi, M., Colvin, G., Liu, Q., Weier, H-U., Johnson, K. and Quesenberry, P.J. (2008) *Conversion potential of marrow cells into lung cells fluctuates with cytokine-induced cell cycle* Stem Cells Devel., **17**, 207-219
39. Dooner, G.J., Colvin, G.A., Dooner, M.S., Johnson, K.W. and Quesenberry, P.J. (2008) *Gene expression fluctuations in murine hematopoietic stem cells with cell cycle progression* J. Cell. Physiol., **214**, 786-795
40. Aliotta, J.M., Pereira, M., Johnson, K.W., de Paza, N., Dooner, M.S., Puente, N., Ayala, C., Brilliant, K., Berza, D., Lee, D., Ramratnam, B., McMillan, P.N., Hixson, D.C., Josic, D. and Quesenberry, P.J. (2010) *Microvesicle entry into marrow cells mediates tissue-specific changes in mRNA by direct delivery of mRNA and induction of transcription* Exp. Hematol., **38**, 233-245
41. Quesenberry, P.J., Dooner, G.J., Del Tatto, M., Colvin, G.A., Johnson, K. and Dooner, M.S. (2010) *Expression of cell cycle-related genes with cytokine-induced cell cycle progression of primitive hematopoietic stem cells* Stem Cells Devel., **19**, 453-460
42. Hu, M.S., Huang, K-J., Li, S., Wu, J-C., Lo, D.D., Hyun, J.S., Chung, M.T., Hu, M., Longaker, M.T. and Lorenz, P. (2014) *Blood-derived mesenchymal stem cells heal calvarial defects and promote wound healing* J. Am. Coll. Surg., **219**, S85-S86
43. Del Carmen Rodriguez, M., Bernad, A. and Aracil, M. (2004) *Interleukin-6 deficiency affects bone marrow stromal precursors, resulting in defective hematopoietic support* Blood, **103**, 3349-3354
44. Franco, S., van de Vrugt, H.J., Fernandez, P., Aracil, M., Arwert, F. and Blasco, M.A. (2004) *Telomere dynamics in Fancg-deficient mouse and human cells* Blood, **104**, 3927-3935
45. Brewer, G. J., Espinosa, J., McIlhane, M. P., Pencek, T. P., Kessler, J. P., Cotman, C., Viel, J. and McManus, D. C. (2001) *Culture and regeneration of human neurons after brain surgery* J. Neurosci. Meth., **107**, 15-23
46. Tatebayashi, Y., Iqbal, K. and Grundke-Iqbal, I. (1999) *Dynamic regulation of expression and phosphorylation of Tau by fibroblast growth factor-2 in neural progenitor cells from adult rat hippocampus* J. Neurosci., **19**, 5245-5254
47. Tatebayashi, Y., Haque, N., Tung, Y-C., Iqbal, K. and Grundke-Iqbal, I. (2004) *Role of tau phosphorylation by glycogen synthase kinase-3 β in the regulation of organelle transport* J. Cell Sci., **117**, 1653-1663
48. Chen, H., Tung, Y-C., Li, B., Iqbal, K. and Grundke-Iqbal, I. (2007) *Trophic factors counteract elevated FGF-2-induced inhibition of adult neurogenesis* Neurobiol. Aging, **28**, 1148-1162
49. Staquicini, F.I., Dias-Neto, E., Li, J., Snyder, E.Y., Sidman, R.L., Pasqualini, R. and Arap, W. (2009) *Discovery of a functional protein complex of netrin-4, laminin γ 1 chain, and integrin α 6 β 1 in mouse neural stem cells* Proc. Natl. Acad. Sci. USA, **106**, 2903-2908
50. Chohan, M.O., Li, B., Blanchard, J., Tung, Y-C., Heaney, A.T., Rabe, A., Iqbal, K. and Grundke-Iqbal, I. (2011) *Enhancement of dentate gyrus neurogenesis, dendritic and synaptic plasticity and memory by a neurotrophic peptide* Neurobiol. Aging, **32**, 1420-1434
51. Abbosh, C., Lawkowski, A., Zaben, M. and Gray, W. (2011) *GalR2/3 mediates proliferative and trophic effects of galanin on postnatal hippocampal precursors* J. Neurochem. (2011) **117**, 425-436
52. He, P. and Shen, Y. (2009) *Interruption of β -catenin signaling reduces neurogenesis in Alzheimer's disease* J. Neurosci., **29**, 6545-6557
53. Nunan, R., Sivasathiaselan, H., Khan, D., Zaben, M. and Gray, W. (2014) *Microglial VPAC1R mediates a novel mechanism of neuroimmune-modulation of hippocampal precursor cells via IL-4 release* Glia, **62**, 1313-1327
54. Yovchev, M.I., Dabeva, M.D. and Oertel, M. (2013) *Isolation, characterization, and transplantation of adult liver progenitor cells* In Methods Mol. Biol., **976**, Stem Cells and Aging: Methods and Protocols, (ed. Tursen, K.) Springer Science+Business Media, LLC pp 37-51

Application Sheet C23; 7th edition, January 2018

Alere Technologies AS

Axis-Shield Density Gradient Media
is a brand of Alere Technologies AS