

# OptiPrep™ Mini-Review MC11

## Isolation of pancreatic stellate cells


- ◆ This Mini-Review primarily provides a complete bibliography of publications that report the use of OptiPrep™ for the purification of pancreatic stellate cells.
- ◆ **OptiPrep™ Application Sheet (C33)** describes the methods available for the purification these cells. It may be accessed from the OptiPrep™ Applications flash-drive or from the following website: [www.axis-shield-density-gradient-media.com](http://www.axis-shield-density-gradient-media.com)
- ◆ This Mini-Review provides short methodological review with a full bibliography.

Pancreatic tissue has been obtained principally from rats, but also from mice and human sources. In the latter case tissue was obtained from patients undergoing resection for malignancy or during the course of surgical treatment for chronic pancreatitis [1]. In its simplest form the disaggregation procedure involves incubation of the finely minced tissue with collagenase, sometimes with pronase, sometimes additionally DNase I is included or added subsequently. With rat tissue, a solution of collagenase (1 mg/5 ml) may first be infused into the pancreas *in situ* via the biliopancreatic duct or thoracic aorta. Subsequently the tissue may be removed and shaken at 37°C for 15 min, prior to further treatment. With or without the infusion, the excised pancreas is often finely minced with scissors and then treated, with collagenase P; for mouse concentrations of 0.03-0.05%, for rat 0.025-0.05%. Concentrations of pronase, if used, are 0.02-0.05% for rat and 0.1% for mouse. DNase I concentrations vary from 0.025% to 0.1% which may be included in the incubation with collagenase or added after a filtration step. Filtration is normally carried out through a nylon mesh. The enzymic digestions are regularly executed at 37°C. For more information on these procedures see refs 1-14.

Occasionally a more elaborate pre-gradient treatment is employed. Kaku et al [15] digested the washed pancreas with 0.12% collagenase for 20 min with shaking; subsequently the digested pancreas was suspended in 10 ml of medium supplemented with 1% BSA and 2 mM EDTA and agitated for 3 min at 37°C; this was repeated three times. A second digestion was carried out with collagenase,  $\alpha$ -chymotrypsin, hyaluronidase and DNase I.

### 2. Density gradient

Occasionally the cells are simply layered on top of a density barrier and the cells recovered from the interface [1-3]. The most commonly used method however is to suspend the cells in culture medium or balanced salt solution (e.g. HBSS) containing 12% (v/v) OptiPrep™; a small volume of medium or HBSS layered on top and then centrifuged at approx. 1,000-1,400 g for 10-30 min [4-13]; stellate cells are recovered from the interface. Sometimes they are suspended in 15% or 17% OptiPrep™ and overlaid with 11.5 or 12.0% OptiPrep™ [14-15], see Figure 1. In the simpler of the two strategies, i.e. suspending the cells in 12% (v/v) OptiPrep™ beneath HBSS (upper part), the stellate cells that float to the interface are essentially adjacent to the original cell suspension layer. This will contain residual disaggregating enzymes, partially broken cells and other pancreatic cells, which do not sediment. In the alternative strategy in which the cells are adjusted to 17% OptiPrep™ and 12% OptiPrep™ layered on top, the stellate cells are separated from the load zone by a clean layer of 12% OptiPrep™ (lower part of Figure 1). It is thus easier to aspirate the stellate cells free from extraneous material by using this two-layer gradient format.



**Figure 1:** Two flotation strategies for purification of pancreatic stellate cells (for details see text).

### 3. References

1. Bachem, M.G., Schneider, E., Groß, H., Weidenbach, H., Schmid, R.M., Menke, A., Siech, M., Beger, H., Grunert, A. and Adler, G. (1998) *Identification, culture, and characterization of pancreatic stellate cells in rat and humans* Gastroenterology, **115**, 421-432
2. Schmid-Kotsas, A., Gross, H-J., Menke, A., Weidenbach, H., Adler, G., Siech, M., Beger, H., Grunert, A. and Bachem, M.G. (1999) *Lipopolysaccharide-activated macrophages stimulate the synthesis of collagen type I and C-fibronectin in cultured pancreatic stellate cells* Am. J. Pathol., **155**, 1749-1758
3. Siech, M., Zhou, Z., Zhou, S., Bair, B., Alt, A., Hamm, S., Gross, H., Mayer, J., Beger, H.G., Tian, X., Kornmann, M. and Bachem, M.G. (2009) *Stimulation of stellate cells by injured acinar cells: a model of acute pancreatitis induced by alcohol and fat (VLDL)* Am. J. Physiol. Gastrointest. Liver Physiol., **297**, G1163-G1171
4. Shek, F.W-T., Benyon, R.C., Walker, F.M., McCrudden, P.R., Pender, S.L.F., Williams, E.J., Johnson, P.A., Johnson, C.D., Bateman, A.C., Fine, D.R. and Iredale, J.P. (2002) *Expression of transforming growth factor- $\beta$ 1 by pancreatic stellate cells and its implications for matrix secretion and turnover in chronic pancreatitis* Am. J. Pathol., **160**, 1787-1798
5. Shek, F.W., Walker, F.M.J., Bateman, A.C., Benyon, R.C., Iredale, J.P. and Fine, D.R. (2003) *A reliable technique for rodent pancreatic stellate cells isolation* Pancreatology, **3**, 209-269
6. Gao, R. and Brigstock, D.R. (2005) *Connective tissue growth factor (CCN2) in rat pancreatic stellate cell function: integrin  $\alpha_5\beta_1$  as a novel CCN2 receptor* Gastroenterology, **129**, 1019-1030
7. Rickmann, M., Vaquero, E.C., Malagelada, J.R. and Molero, X. (2007) *Tocotrienols induce apoptosis and autophagy in rat pancreatic stellate cells through the mitochondrial death pathway* Gastroenterology, **132**, 2518-2532
8. González, A.M., Garcia, T., Samper, E., Rickmann, M., Vaquero, E.C. and Molero, X. (2011) *Assessment of the protective effects of oral tocotrienols in arginine chronic-like pancreatitis* Am. J. Physiol. Gastrointest. Liver Physiol., **301**, G846-G855
9. Armstrong, T., Packham, G., Murphy, L.B., Bateman, A.C., Conti, J.A., Fine, D.R., Johnson, C.D., Benyon, R.C. and Iredale, J.P. (2004) *Type I collagen promotes the malignant phenotype of pancreatic ductal adenocarcinoma* Clin. Cancer Res., **10**, 7427-7437
10. José, A., Rovira-Rigau, M., Luna, J., Giménez-Alejandre, M., Vaquero, E., García de la Torre, B., Andreu, D., Alemany, R. and Fillat, C. (2014) *A genetic fiber modification to achieve matrix-metalloprotease-activated infectivity of oncolytic adenovirus* Journal of Control. Release, **192**, 148-156
11. Ulmasov, B., Xu, Z., Tetri, L.H., Inagami, T. and Neuschwander-Tetri, B.A. (2009) *Protective role of angiotensin II type 2 receptor signaling in a mouse model of pancreatic fibrosis* Am. J. Physiol. Gastrointest. Liver Physiol., **296**, G284-G294
12. Lawencia, C., Charrier, A., Huang, G. and Brigstock, D.R. (2009) *Ethanol-mediated expression of connective tissue growth factor (CCN2) in mouse pancreatic stellate cells* Growth Factors, **27**, 91-99
13. Charrier, A., Chen, R., Chen, L., Kemper, S., Hattori, T., Takigawa, M. and Brigstock, D.R. (2014) *Connective tissue growth factor (CCN2) and microRNA-21 are components of a positive feedback loop in pancreatic stellate cells (PSC) during chronic pancreatitis and are exported in PSC-derived exosomes* J. Cell Commun. Signal., **8**, 147-156
14. Kaku, T., Oono, T., Zhao, H., Gibo, J., Kawabe, K., Ito, T. and Takayanagi, R. (2007) *IS-741 attenuates local migration of monocytes and subsequent pancreatic fibrosis in experimental chronic pancreatitis induced by dibutyltin dichloride in rats* Pancreas, **34**, 299-309
15. Yu, F-X., Su, L-F., Dai, C-L., Wang, Y., Teng Y-Y., Fu, J-H., Zhang, Q-Y. and Tang, Y-H. (2015) *Inhibition of pancreatic stellate cell activity by adipose-derived stem cells* Hepatobiliary Pancreat. Dis. Int., **14**, 215-221
16. Erkan, M., Adler, G., Apte, M.V., Bachem, M.G., Buchholz, M., Detlefsen, S., Esposito, I., Friess, H., Gress, T.M., Habisch, H-J. et al (2012) *StellaTUM: current consensus and discussion on pancreatic stellate cell research* Gut, **61**, 172-178

Mini-Review MC11: 3<sup>rd</sup> edition, January 2017

Alere Technologies AS

Axis-Shield Density Gradient Media  
is a brand of Alere Technologies AS