

OptiPrep™ Application Sheet S49

Isolation of optical disks from retinal rod outer segments

- ◆ OptiPrep™ is a 60% (w/v) solution of iodixanol in water, density = 1.32 g/ml
- ◆ To access other Application Sheets referred to in the text return to the Subcellular Membranes Index; key Ctrl “F” and type the S-Number in the Find Box.
- ◆ Important technical notes, information regarding alternative methodologies and membrane analysis are contained in the “Technical Notes and Review” section (Section 5).
- ◆ Section 5.3 includes information on some variations to the given format and some of the more recent references

1. Background

The use of OptiPrep™ for the isolation of mouse retinal rod outer segments (ROSs) was first published in 1998 [1]. After removal of the retinas, they are repeatedly vortexed in a Ringer’s solution containing iodixanol (usually, but not always, 8% v/v/ OptiPrep™) and centrifuged at low speed. During this procedure the ROS s are released into the supernatants and subsequently banded in a two-step 10%, 18% (v/v) OptiPrep™ gradient. The method was extended without modification to the isolation of ROSs from human tissue [2]. Howes et al [3] and Calvert et al [4] replaced the two-layer iodixanol gradient with a three layer one for the purification of the ROSs, although the densities of the three layers were slightly different and Ringer’s solution was replaced by Locke’s solution.

More recently Liang et al [5] introduced a continuous iodixanol gradient for the purification of the ROSs and subsequently, after osmotic lysis of the ROSs, the optical disks were recovered on a second continuous iodixanol gradient. For convenience, the methodology described in this Application Sheet is adapted from ref 5 but the extraction and gradient variations for ROS purification are presented in Section 5.3.

2. Solutions required (see Section 5.1)

- A. OptiPrep™
- B. Ringers Solution: 130 mM NaCl, 3.6 mM KCl, 2.4 mM MgCl₂, 1.2 mM CaCl₂, 0.02 mM EDTA, 10 mM Hepes-NaOH, pH 7.4
- C. Lysis buffer: 2 mM Tris-HCl, pH 7.4

3. Centrifuge requirements

Swinging-bucket rotor for a high-speed centrifuge with a tube volume of approx 15 ml

4. Protocol

1. Make an 8% (v/v) OptiPrep™ solution by diluting 8 ml of Solution A to 100 ml with Solution B.
2. Suspend 12 mouse retinas in 0.12 ml of this solution. Carry out this and the following steps at room temperature.
3. Vortex for 1 min and then centrifuge at 200 g for 1 min.
4. Carefully recover the supernatant (without disturbing the pellet) and retain.
5. Resuspend the pellet in 0.12 ml of the 8% OptiPrep™ solution and repeat steps 3-4.
6. Repeat this procedure five or six times.

Keep the following stock solutions at 4°C:
 500 mM Hepes (free acid): 11.9 g per 100 ml water
 10 mM EDTA (Na₂•2H₂O): 0.372 g per 100 ml water
 1 M NaCl: 5.84 g per 100 ml of water
 100 mM KCl: 0.745 g per 100 ml of water
 100 mM MgCl₂•6H₂O: 2.03 g per 100 ml of water
 100 mM CaCl₂•2H₂O: 1.47 g per 100 ml water
 100 mM Tris (free base): 1.21 g per 100 ml

Solution B: Mix together the following volumes of the above stocks, **except the Tris** (in order), 4 ml, 0.4 ml, 26 ml, 7.2 ml, 4.8 ml and 2.4 ml; adjust to pH 7.4 with NaOH and make up to 200 ml

Solution C: Add 4 ml of Tris stock solution to 100 ml of water; adjust to pH 7.4 with 1 M HCl; make up to 200 ml

7. Prepare two gradient solutions of 10% and 30% (v/v) OptiPrep™ by diluting 1 ml and 3 ml of Solution A to 10 ml with Solution B respectively.
8. Using a two-chamber gradient maker or a Gradient Master™ prepare a 10-12 ml continuous gradient from these two solutions in tubes for the swinging-bucket rotor (see [Section 5.2](#)).
9. Layer the combined supernatants (approx. 1.5 ml) on top and centrifuge at 26,500 g for 30 min.
10. Collect the band of ROSs (about two thirds of the way from the top of the gradient).
11. Dilute with 3 vol. of Solution B and centrifuge at 500 g for 3 min to pellet the nuclei.
12. Harvest the ROSs by centrifuging at 26,500 g for 30 min.
13. To prepare the optical disks, lyse the ROSs by suspending them in 2 ml of Solution C and maintain them at 0°C for 15 h.
14. Prepare two gradient solutions of 15% and 40% (v/v) OptiPrep™ from 1.5 ml and 4 ml of Solution A, each diluted to 10 ml with Solution B respectively.
15. Using a two-chamber gradient maker or a Gradient Master™ prepare a 10-12 ml continuous gradient from these two solutions in tubes for the swinging-bucket rotor (see [Section 5.2](#)).
16. Layer the lysed ROSs on top of the gradient and centrifuge at 26,500 g for 30 min.
17. Harvest the disks that band about two thirds of the way from the top of the gradient

5. Technical Notes and Review

5.1 Homogenization media and gradient solutions

All iodixanol solutions in this Optiprep™ Application Sheet are prepared simply by diluting Optiprep™ with Ringers Solution, which will mean that the levels of KCl, MgCl₂, CaCl₂, EDTA and buffer will be reduced. If it is considered important that these levels should be maintained in all density solutions then a working solution of 50% (w/v) iodixanol should be first prepared by diluting 5 vol. of Optiprep™ with 1 vol. of 130 mM NaCl, 21.6 mM KCl, 14.4 mM MgCl₂, 7.2 mM CaCl₂, 0.12 mM EDTA, 60 mM Hepes-NaOH, pH 7.4. The working solution would thus contain the same concentrations of KCl, MgCl₂, CaCl₂, EDTA and Hepes as the Ringers buffer. Further dilutions with Ringers buffer should then need to be adjusted appropriately from those given in the protocol to give gradient solutions of the correct density. Note that if the concentration of NaCl in the diluent were also raised then the osmolality of the 50% (w/v) iodixanol will be unacceptably high. Strategies for preparing working solutions for mammalian tissues and cells are given in [Application Sheet S01](#).

5.2 Construction of gradients

If neither of these gradient-making devices is available then a continuous gradient can be prepared by diffusion of a discontinuous gradient. For example, the 10-30% (v/v) Optiprep™ gradient (step 7) might be generated from a discontinuous gradient of equal volumes of 10%, 18%, 24% and 30% Optiprep™. Likewise the 15-40% (v/v) Optiprep™ gradient (step 14) might be generated from a discontinuous gradient of equal volumes of 15%, 21%, 27%, 33% and 40% Optiprep™. For more information on gradient construction see [Application Sheet S03](#).

5.3 Variations in gradient and centrifugation format

5.3.1 Continuous gradients

Continuous iodixanol gradients normally conform to the above protocol [5-13], although Peshenko et al [14,15] used a 20-30% (v/v) OptiPrep™ gradient and a much higher *g*-force of 75,000 *g*, for 40 min.

5.3.2 Discontinuous gradients

The commonly used discontinuous gradient format involves layering the extracted ROSs in 8% (v/v) Optiprep™ over a two layers of 10% and 18% (v/v) Optiprep™; the centrifugation conditions are however somewhat diverse: 3,300 *g* for 10 min [1,2], 50,000 *g* for 1 h [16] 70,000 *g* for 1 h [17]. Variants in the density of the gradient layers include 8.3% and 20% with the sample in 6.7%, centrifuged at only 1,425 *g*, for 15 min [4]. Gradients comprising three layers: 10, 20 and 30% (sample in 8%), after centrifugation at 17,000 *g* for 50 min the ROSs banded at the 20%/30% interface [3] and four layers: 8%, 12%, 16% and 20% (sample in 2%) [18,19].

- ◆ Occasionally after repeated extractions with 8% (v/v) Optiprep™ no gradient is used [20].
- ◆ Refs 21-23 are meetings abstracts, reporting the use of Optiprep™ but no details are provided.
- ◆ **Refs 24-32 all report the use of iodixanol gradients prepared from OptiPrep™ for ROS isolation using methods identical or very similar to those described above**

6. References

1. Tsang, S.H., Burns, M.E., Calvert, P.D., Gouras, P., Baylor, D.A., Goff, S.P. and Arshavsky, V.Y. (1998) *Role for the target enzyme in deactivation of photoreceptor G protein in vivo* Science, **282**, 117-121
2. Mata, N.L., Weng, J. and Travis, G.H. (2000) *Biosynthesis of a major lipofuscin fluorophore in mice and humans with ABCR-mediated retinal and macular degeneration* Proc. Natl. Acad. Sci. USA, **97**, 7154-7159
3. Howes, K.A., Pennesi, M.E., Sokal, I., Church-Kopish, J., Schmidt, B., Margolis, D., Frederick, J.M., Rieke, F., Palczewski, K., Wu, S.M., Detwiler, P.B. and Baehr, P. (2002) *GCAP1 and rescues rod photoreceptor response in GCAP1/GCAP2 knockout mice* The EMBO J., **21**, 1545-1554
4. Calvert, P.D., Govardovskii, V.I., Krasnoperova, N., Anderson, R.E., Lem, J. and Makino, C.L. (2001) *Membrane protein diffusion sets the speed of rod phototransduction* Nature, **411**, 90-94
5. Liang, Y., Fotiadis, D., Filipek, S., Saperstein, D. A., Palczewski, K. and Engel, A. (2003) *Organization of the G protein-coupled receptors rhodopsin and opsin in native membranes* J. Biol. Chem., **278**, 21655-21662
6. Fotiadis, D., Liang, Y., Filipek, S., Saperstein, D. A., Engel, A. and Palczewski, K. (2004) *The G protein-coupled receptor rhodopsin in the native membrane* FEBS Lett., **564**, 281-288
7. Liang, Y., Fotiadis, D., Filipek, S., Saperstein, D.A., Engel, A. and Palczewski, K. (2004) *Rhodopsin signaling and organization in heterozygote rhodopsin knockout mice* J. Biol. Chem., **279**, 48189-48196
8. Maeda, A., Maeda, T., Imanishi, Y., Kuksa, V., Alekseev, A., Bronson, J.D., Zhang, H., Zhu, L., Sun, W., Saperstein, D.A., Ricke, F., Baehr, W. and Palczewski, K. (2005) *Role of photoreceptor-specific retinal dehydrogenase in the retinoid cycle in vivo* J. Biol. Chem., **280**, 18822-18832
9. Saperstein, D.A., Fotiadis, D., Liang, Y., Filipek, S., Palczewski, K. and Engel, A. (2003) *The structure of murine outer segment disk membranes using atomic force microscopy* Invest. Ophthalmol Vis. Sci., **44**, E-abstract, 3175
10. Nickell, S., Park, P.S-H., Baumeister, W. and Palczewski, K. (2007) *Three-dimensional architecture of murine outer rod segments determined by cryoelectron tomography* J. Cell Biol., **177**, 917-925
11. Mustafi, D., Kevany, B.M., Genoud, C., Okano, K., Cideciyan, A.V., Sumaroka, A., Roman, A.J., Jacobson, S.G., Engel, A., Adams, M.D. and Palczewski, K. (2011) *Defective photoreceptor phagocytosis in a mouse model of enhanced S-cone syndrome causes progressive retinal degeneration* FASEB J. **25**, 3157-3176
12. Ziccardi, L., Vijayarathy, C., Bush, R.A. and Sieving, P.A. (2012) *Loss of retinoschisin (RS1) cell surface protein in maturing mouse rod photoreceptors elevates the luminance threshold for light-driven translocation of transducin but not arrestin* J. Neurosci., **32**, 13010-13021
13. Gilliam, J.C., Chang, J.T., Sandoval, I.M., Zhang, Y., Li, T., Pittler, S.J., Chiu, W. and Wensel, T.G. (2012) *Three-dimensional architecture of the rod sensory cilium and its disruption in retinal neurodegeneration* Cell. **151**, 1029-1041
14. Peshenko, I.V., Olshevskaya, E.V., Savchenko, A.B., Karan, S., Palczewski, K., Baehr, W. and Dizhoor, A.M. (2011) *Enzymatic properties and regulation of the native isozymes of retinal membrane guanylyl cyclase (RetGC) from mouse photoreceptors* Biochemistry, **50**, 5590-5600

15. Peshenko, I.V., Olshevskaya, E.V., Azadi, S., Molday, L.L., Molday, R.S. and Dizhoor, A.M. (2011) *Retinal degeneration 3 (RD3) protein inhibits catalytic activity of retinal membrane guanylyl cyclase (RetGC) and its stimulation by activating proteins* Biochemistry, **50**, 9511– 9519
16. Dizhoor, A.M. Woodruff, M.L., Olshevskaya, E.V., Cilluffo, M.C., Cornwall, M.C., Sieving, P.A. and Fain, G.L. (2008) *Night blindness and the mechanism of constitutive signaling of mutant G90D rhodopsin* J. Neurosci., **28**, 11662-11672
17. Burns, M.E., Mendez, A., Chen, C-K., Almuete, A., Quillinan, N., Simon, M.I., Baylor, D.A. and Chen, J. (2006) *Deactivation of phosphorylated and nonphosphorylated rhodopsin by arrestin splice variants* J. Neurosci., **26**, 1036-1044
18. Krispel, C.M., Chen, D., Melling, N., Chen, Y-J., Martemyanov, K.A., Quillinan, N., Arshavsky, V.Y., Wensel, T.G., Chen, C-K. and Burns, M.E. (2006) *RGS expression rate-limits recovery of rod photoresponses* Neuron, **51**, 409-416
19. Gross, O.P. and Burns, M.E. (2010) *Control of rhodopsin's active lifetime by arrestin-1 expression in mammalian rods* J. Neurosci., **30**, 3450 –3457
20. Nair, K.S., Hanson, S.M., Mendez, A., Gurevich, E.V., Kennedy, M.J., Shestopalov, V.I., Vishnivetskiy, S.A., Chen, J., Hurley, J.B., Gurevich, V.V. and Slepak, V.Z. (2005) *Light-dependent redistribution of arrestin in vertebrate rods is an energy-independent process governed by protein-protein interactions* Neuron, **46**, 555-567
21. Wang, Q., Hu, G., Leitges, M. and Wensel, T.G. (2006) *Phosphorylation of phototransduction GAP RGS9-1 depends on the α isoform of protein kinase C* Invest. Ophthalmol. Vis. Sci., **47**, E-Abstr 822
22. Finnemann, S.C. and Chang, Y. (2006) *Regulation of RPE phagocytosis by integrin receptor-tetraspanin surface membrane domains* Invest. Ophthalmol. Vis. Sci., **47**, E-Abstr 822
23. Coleman, J.A., Djajadi, H.R., Molday, L.L. and Molday, R.S. (2012) *Disruption of the P₄-ATPase aminophospholipid flippase Atp8a2 gene suggests a role for phosphatidylserine in photoreceptor outer segments* Invest. Ophthalmol. Vis. Sci., **53**, Abstr. 743- 620
24. Moaven, H., Koike, Y., Jao, C.C., Gurevich, V.V., Langen, R. and Chen, J. (2013) *Visual arrestin interaction with clathrin adaptor AP-2 regulates photoreceptor survival in the vertebrate retina* Proc. Natl. Acad. Sci. USA, **110**, 9463–9468
25. Skiba, N.P., Spencer, W.J., Salinas, R.Y., Lieu, E.C., Thompson, J.W. and Arshavsky, V.Y. (2013) *Proteomic identification of unique photoreceptor disc components reveals the presence of PRCD, a protein linked to retinal degeneration* J. Proteome Res., **12**, 3010-3018
26. Coleman, J.A., Zhu, X., Djajadi, H.R., Molday, L.L., Smith, R.S., Libby, R.T., John, S.W.M. and Molday, R.S. (2014) *Phospholipid flippase ATP8A2 is required for normal visual and auditory function and photoreceptor and spiral ganglion cell survival* J. Cell Sci., **127**, 1138–1149
27. Wensel, T.G. and Gilliam, J.C. (2015) *Three-dimensional architecture of murine rod cilium revealed by cryo-EM* In Methods Mol. Biol., **1271**, Rhodopsin: Methods and Protocols (ed. Jastrzebska, B.), Springer Science+Business Media, New York, pp 267-292
28. Rakshit, T., Senapati, S., Sinha, S., Whited, A.M. and Park, P.S-H. (2015) *Rhodopsin forms nanodomains in rod outer segment disc membranes of the cold-blooded Xenopus laevis* PLoS One, **10**: e0141114
29. Calzia, D., Panfoli, I., Heinig, N., Schumann, U., Ader, M., Traverso, C.E., Funk, R.H.W. and Roehlecke, C. (2016) *Impairment of extramitochondrial oxidative phosphorylation in mouse rod outer segments by blue light irradiation* Biochimie, **125**, 171-178
30. He, F., Agosto, M.A., Anastassov, I.A., Tse, D.Y., Wu, S.M. and Wensel, T.G. (2016) *Phosphatidylinositol-3-phosphate is light-regulated and essential for survival in retinal rods* Sci. Rep., **6**: 26978
31. Berry, J., Frederiksen, R., Yao, Y., Nymark, S., Chen, J. and Cornwall, C. (2016) *Effect of rhodopsin phosphorylation on dark adaptation in mouse rods* J. Neurosci., **36**, 6973– 6987
32. Peshenko, I.V., Olshevskaya, E.V. and Dizhoor, A.M. (2016) *Functional study and mapping sites for interaction with the target enzyme in retinal degeneration 3 (RD3) protein* J. Biol. Chem., **291**, 19713–19723