

OptiPrep™ Application Sheet S55

Isolation of peroxisomes from yeast spheroplasts

- ◆ OptiPrep™ is a 60% (w/v) solution of iodixanol in water, density = 1.32 g/ml
- ◆ To access other Application Sheets referred to in the text return to the Subcellular Membranes Index; key Ctrl “F” and type the S-Number in the Find Box.
- ◆ Recent published papers reporting the use of iodixanol gradients are given in Section 7.

1. Background

Wild-type yeast strains contain rather few peroxisomes, however they can be induced to proliferate and there are a number of peroxisome gene-deletion strains of the organism [1,2]. Yeast systems are therefore a popular means of investigating peroxisome development.

Although Percoll® has been used for peroxisome isolation from mammalian cells and yeast spheroplasts, the organelles are contaminated by both Golgi and smooth ER membranes, which have similar low densities in this medium. In both iodixanol and Nycodenz® however peroxisomes are denser than any other organelle in the light mitochondrial fraction and contamination by ER and Golgi is never a problem. Nycodenz® gradients have been widely used for the isolation of peroxisomes, but in the last eight years there has been increasing use of OptiPrep™ for both mammalian and yeast peroxisomes. Although Nycodenz® has remained a popular choice for yeast peroxisomes, the availability of iodixanol as a 60% (w/v) solution (OptiPrep™) makes gradient solution preparation much easier than with Nycodenz®. Iodixanol and Nycodenz® solutions of the same % (w/v) concentration have almost identical densities, but only iodixanol permits the banding of peroxisomes under more or less isoosmotic conditions and it is only the OptiPrep™ option that is presented in this Application Sheet.

Spheroplasts are homogenized in a sorbitol-containing medium using a Dounce homogenizer. Some variation exists in the sorbitol concentration in this medium; Watkins et al [2] used 0.6 M, while Crane et al [3] used 1.0 M. The iodixanol methodology is adapted from refs 4-5, 7-10.

2. Solutions required (see Box)

- Spheroplast wash solution: 1.2 M sorbitol, 20 mM phosphate buffer, pH 7.4
- Spheroplast lysis buffer: 0.6 M sorbitol, 1 mM EDTA, 1 mM KCl, 0.1% (v/v) ethanol, 5 mM MES-NaOH pH 6.0
- OptiPrep™
- OptiPrep™ diluent: 18% (w/v) sucrose, 3 mM EDTA, 3 mM KCl, 0.3% ethanol, 15 mM MES-NaOH, pH 6.0
- Gradient diluent: 18% (w/v) sucrose, 1 mM EDTA, 1 mM KCl, 0.1% ethanol, 5 mM MES-NaOH, pH 6.0

Add protease inhibitors to solutions as required.

3 Preparation of gradient solutions

Mix 2 vol. of OptiPrep™ with 1 vol. of Solution D to produce a 40% (w/v) iodixanol stock solutions and then dilute further with Solution E to produce solutions of 2.25% and 24% (w/v) iodixanol (see Section 6, Notes 1 and 2)

Keep the following stock solutions at 4°C:

1 M KCl	7.45 g per 100 ml water
1 M MES, purchased ready-made solution	
100 mM KH ₂ PO ₄	1.36 g per 100 ml water
100 mM K ₂ HPO ₄ •3H ₂ O	2.28 g per 100 ml water
100 mM EDTA (Na ₂ •2H ₂ O)	3.72 g per 100 ml water

Solution A: Mix 40 ml of K₂HPO₄ stock with 10 ml of KH₂PO₄ stock; check the pH is 7.4 and adjust if necessary. Dissolve 21.8 g of sorbitol in 50 ml of water, add 20 ml of the phosphate buffer and make up to 100ml

Solution B: Dissolve 10.9 g of sorbitol in 50 ml of water; add 0.1 ml, 1 ml and 0.5 ml of the KCl, EDTA and MES stocks respectively; adjust to pH 6.0 with 0.1 M NaOH; add 0.1 ml of ethanol and make up to 100 ml.

Solution D: Dissolve 18 g of sucrose in 50 ml of water; add 0.3 ml, 3 ml and 1.5 ml of the KCl, EDTA and MES stocks respectively; adjust to pH 6.0 with 0.1 M NaOH; add 0.3 ml of ethanol and make up to 100 ml.

Solution E: As Solution C using one third of the volumes of KCl, EDTA, MES and alcohol

4 Ultracentrifuge rotor requirements

Vertical rotor with a tube capacity of approx 36 ml tubes (e.g. Beckman VTi50) with appropriate sealed tubes (see Section 6, Note 3).

5. Protocol

Carry out all operations at 0-4°C

1. Prepare spheroplasts from 1 litre of yeast culture, grown in YPD medium ($OD_{600} = 0.5-1.0$).
2. Suspend the spheroplasts in Solution A (10-15 ml).
3. Centrifuge in an 8x50 ml fixed-angle rotor (high-speed centrifuge) at 4000 g for 5 min.
4. Remove the supernatant and repeat steps 1 and 2.
5. Suspend the spheroplasts in 35 ml of Solution B and homogenize in a tight-fitting Dounce homogenizer (Wheaton Type A), using 10 up-and-down strokes of the pestle.
6. Centrifuge the homogenate at 1,500 g for 10 min.
7. Using a syringe and metal cannula aspirate and retain the supernatant on ice (see Section 6, Note 4).
8. Resuspend the pellet in 35 ml of Solution B and repeat steps 5-7.
9. Combine the two supernatants and centrifuge at 25,000 g for 30 min.
10. Resuspend the light mitochondrial pellet in 6 ml of Solution B using 10 gentle strokes of the pestle of the small volume loose-fitting (Wheaton Type B) Dounce homogenizer (see Section 6, Notes 5 and 6).
11. Using a two chamber gradient maker or Gradient Master™ prepare 30 ml linear gradients from equal volumes of the 2.25 and 24% (w/v) iodixanol in tubes for the vertical rotor (see Section 6, Notes 7 and 8).
12. Underlayer the gradient with 1.0 ml of 0.5 ml of OptiPrep™.
13. Layer the sample on top of the gradient to fill the tube and seal it (see Section 6, Note 9).
14. Centrifuge at 40,000 g for 90 min. Use a controlled acceleration and deceleration programs to ensure a smooth reorientation of the gradient. If these are not available, turn off the brake below 2000 rpm.
15. Collect the gradient dense end first in 0.5 ml fractions; the peroxisomes band close to the bottom of the gradient. For more information on gradient harvesting and analysis see respectively Application Sheets S08 and S09.

6. Notes

1. In the gradient EDTA was omitted in the method of ref 7.
2. Originally the gradient covered a wider density range 15-36% (w/v) iodixanol [4,5], but more recently the 2.25-24% or 2.25-22.5% [7-14] iodixanol gradients have become more popular. Discontinuous iodixanol gradients are less often [15,16]; the light mitochondrial fraction was adjusted to 23.5% (w/v) iodixanol and layered over 35% (w/v) iodixanol and centrifuged at 110,000 g for 2 h. The peroxisomes band at the interface.
3. The use of a vertical rotor is very common in many organelle purifications. The short sedimentation path length of the rotor means that the particles reach their banding density very

quickly and the low hydrostatic pressure in the gradient preserves organelle integrity. A smaller rotor such as the VTi65.1 is permissible. If a vertical rotor is not available, fixed-angle rotor or swinging-bucket rotors may be substituted, but the centrifugation times will need to be increased.

4. Metal cannulas (i.d. approx 0.8 mm) can be obtained from most surgical instrument supply companies.
5. To avoid damage to the delicate organelles only use very gentle strokes of the pestle.
6. A post-nuclear supernatant was used instead of a light mitochondrial fraction according to ref 5.
7. If neither of these devices is available then first construct a discontinuous gradient from equal volumes of 15%, 25%, 30% and 36% (w/v) iodixanol and allow them to diffuse. For more information on preparing gradients **see Application Sheet S03**.
8. For smaller rotors scale down all volumes proportionately.
9. There are a variety of sealed tubes that are commercially available but the easiest sealed tubes to use are Beckman Optiseal™ tubes.

7. Recent publications

In a Cold Spring Harbor Protocols publication Cramer et al [17] provided a detailed and interesting comparison of sucrose gradients, iodixanol gradients in which the OptiPrep™ was diluted with a simple buffer and iodixanol gradients in which the OptiPrep™ was diluted with a standard buffered sucrose homogenization medium and found that the latter clearly gave superior results. Wróblewska et al [18] using an iodixanol gradient containing 18% (w/v) sucrose, based on the methods described by Cramer [17] were able to resolve not only very pure peroxisomes, but also a population of vesicles that contained a subset of peroxisomal membrane proteins. Effelsberg et al [19] used a 15.5-36% (w/v) iodixanol gradient, also containing 18% (w/v) sucrose. Mindthoff et al [20] and Effelsberg et al [21], used the 2.25-25% (w/v) iodixanol gradient previously described by Grunau et al [9].

8. References

1. Thieringer, R., Shio, H., Han, Y., Cohen, G. and Lazarow, P. B. (1991) *Peroxisomes in Saccharomyces cerevisiae: immunofluorescence analysis and import of catalase A into isolated peroxisomes* Mol. Cell. Biol., **11**, 510-522
2. Watkins, P. A., Lu, J-F., Steinberg, S. J., Gould, S. J., Smith, K. D. and Braiterman, L. T. (1998) *Disruption of the Saccharomyces cerevisiae FAT1 gene decreases very long-chain fatty acyl-CoA synthetase activity and elevates intracellular very long-chain fatty acid concentrations* J. Biol. Chem., **273**, 18210-18219
3. Crane, D. I., Kalish, J. E. and Gould, S. J. (1994) *The Pichia pastoris PAS4 gene encodes a ubiquitin-conjugating enzyme required for peroxisome assembly* J. Biol. Chem., **269**, 21835-21844
4. Einwachter, H., Sowinski, S., Kunau, W-H. and Schliebs, W. (2001) *Yarrowia lipolytica Pex20p, Saccharomyces cerevisiae Pex18p/Pex 21p and mammalian Pex5pL fulfil a common function in the early steps of the peroxisomal PTS2 import pathway* EMBO Rep., **2**, 1035-1039
5. Schäfer, A., Kerssen, D., Veenhuis, M., Kunau, W-H. and Schliebs, W. (2004) *Functional similarity between the peroxisomal PTS2 receptor binding protein Pex18p and the N-terminal half of the PTS1 receptor Pex5p* Mol. Cell Biol., **24**, 8895-8906
6. Kerssen, D., Hambruch, E., Klaas, W., Platta, H.W., de Kruijff, B., Erdmann, R., Kunau, W-H. and Schliebs, W. *Membrane association of the cycling peroxisome import receptor Pex5p* J. Biol. Chem., **281**, 27003-27015
7. Platta, H.W., Grunau, S., Rosenkrantz, K., Girzalsky, W. and Erdmann, R. (2005) *Functional role of the AAA peroxins in dislocation of the cycling PTS1 receptor back to the cytosol* Nat. Cell Biol., **7**, 817-822
8. Thoms, S., Debelyy, M.O., Nau, K., Meyer, H.E. and Erdmann, R. (2008) *Lpx1p is a peroxisomal lipase required for normal peroxisome morphology* FEBS J., **275**, 504-514
9. Grunau, S., Mindthoff, S., Rottensteiner, H., Sormunen, R.T., Hiltunen, J.K., Erdmann, R. and Antonenkov, V.D. (2009) *Channel-forming activities of peroxisomal membrane proteins from the yeast Saccharomyces cerevisiae* FEBS J., **276**, 1698-1708
10. Antonenkov, V.D., Mindthoff, S., Grunau, S., Erdmann, R. and Hiltunen, J.K. (2009) *An involvement of yeast peroxisomal channels in transmembrane transfer of glyoxylate cycle intermediates* Int., J. Biochem. Cell Biol., **41**, 2546-2554
11. Welker, S., Rudolph, B., Frenzel, E., Hagn, F., Liebisch, G., Schmitz, G., Scheuring, J., Kerth, A., Blume, A., Weinkauff, S., Haslbeck, M., Kessler, H. and Buchner, J. (2010) *Hsp12 is an intrinsically unstructured stress protein that folds upon membrane association and modulates membrane function* Mol. Cell, **39**, 507-520

12. Grunau, S., Lay, D., Mindthoff, S., Platta, H.W., Girzalsky, W., Just, W.W. and Erdmann, R. (2011) *The phosphoinositide 3-kinase Vps34p is required for pexophagy in Saccharomyces cerevisiae* Biochem. J. **434**, 161–170
13. Debelyy, M.O., Platta, H.W., Saffian, D., Hensel, A., Thoms, S., Meyer, H.E., Warscheid, B., Girzalsky, W. and Erdmann, R. (2011) *Ubp15p, a ubiquitin hydrolase associated with the peroxisomal export machinery* J. Biol. Chem., **286**, 28223–28234
14. Oeljeklaus, S., Reinartz, B.S., Wolf, J., Wiese, S., Tonillo, J., Podwojski, K., Kuhlmann, K., Stephan, C., Meyer, H.E., Schliebs, W., Brocard, C., Erdmann, R. and Warscheid, B. (2012) *Identification of core components and transient interactors of the peroxisomal importomer by dual-track stable isotope labeling with amino acids in cell culture analysis* J. Proteome Res. 2012, 11, 2567–2580
15. Nyathi, Y., De Marcos Lousa, C., van Roermund, C.W., Wanders, R.J.A., Johnson, B., Baldwin, S.A., Theodoulou, F.L. and Baker, A. (2010) *The Arabidopsis peroxisomal ABC transporter, Comatose, complements the Saccharomyces cerevisiae pxa1 pxa2Δ mutant for metabolism of long-chain fatty acids and exhibits fatty acyl-CoA-stimulated ATPase activity* J. Biol. Chem., **285**, 29892–29902
16. Nyathi, Y., Zhang, X., Baldwin, J.M., Bernhardt, K., Johnson, B., Baldwin, S.A., Theodoulou, F.L. and Baker, A. (2012) *Pseudo half-molecules of the ABC transporter, COMATOSE, bind Pex19 and target to peroxisomes independently but are both required for activity* FEBS Lett., **586**, 2280–2286
17. Cramer, J., Effelsberg, D., Girzalsky, W. and Erdmann, R. (2015) *Isolation of peroxisomes from yeast* Cold Spring Harb. Protoc; doi:10.1101/pdb.top074500
18. Wróblewska, J.P., Cruz-Zaragoza, L.D., Yuan, W., Schummer, A., Chuartzman, S.G., de Boer, R., Oeljeklaus, S., Schuldiner, M. et al (2017) *Saccharomyces cerevisiae cells lacking Pex3 contain membrane vesicles that harbor a subset of peroxisomal membrane proteins* BBA – Mol. Cell Res., **1864**, 656–1667
19. Effelsberg, D., Cruz-Zaragoza, L.D., Tonillo, J., Schliebs, W. and Erdmann, R. (2015) *Role of Pex21p for piggyback import of Gpd1p and Pnc1p into peroxisomes of Saccharomyces cerevisiae* J. Biol. Chem., **290**, 25333–25342
20. Mindthoff, S., Grunau, S., Steinfert, L.L., Girzalsky, W., Hiltunen, J.K., Erdmann, R. and Antonenkov, V.D. (2016) *Peroxisomal Pex11 is a pore-forming protein homologous to TRPM channels* Biochim. Biophys. Acta, **1863**, 271–283
21. Effelsberg, D., Cruz-Zaragoza, L.D., Schliebs, W. and Erdmann, R. (2016) *Pex9p is a new yeast peroxisomal import receptor for PTS1-containing proteins* J. Cell Sci., **129**, 4057-4066

Application Sheet S55; 5th edition, February 2018

Alere Technologies AS

Axis-Shield Density Gradient Media
is a brand of Alere Technologies AS