

### **Growing up of cells**

Streak out His-tagged Xis strain from frozen supply

The next day pick a well isolated colony and start a 5 ml overnight per 500 ml to be prepared (2L suggested)

Early the next day, subculture overnight 1:100 into 500 ml LB plus appropriate antibiotic

Grow to mid to late log at 37C OD600=0.5-0.6 (takes about 4 hours)

Induce with a final IPTG concentration of 0.5 mM

Shake at room temperature for four more hours

Spin cells at least at 5000rpm for 20 minutes

Decant all LB

Freeze pellet overnight

### **Cell extracts**

\*all steps performed on ice\*

Thaw pellet on ice for 30-60 minutes

Resuspend pellet with 2ml of extraction buffer

50 mM Tris, pH 7.4, 10% sucrose, 0.8M KCl

Transfer to large homogenizer (25 mL capacity)

Add protease inhibitors

5ul PMSF 50mM in EtOH - made fresh

4ul leupeptin 5mg/ml

2ul pepstatin A 25 mg/ml in EtOH

1ul soybean trypsin inhibitor 100mg/ml

Add 1/20 total volume freshly made lysozyme (10 mg/ml in 0.25M Tris, pH 7.4)

Homogenize until even, watery consistency about 5 – 10 minutes

Incubate on ice 45 minutes

Spin at 45,000 rpm for 45 minutes at 4C

Transfer supernatant to a new eppendorf tube, take an aliquot, and freeze both in an ethanol/dry ice bath or liquid nitrogen

Transfer pellet back into homogenizer and add 0.5 mL extraction buffer – rehomogenize and spin again

Save the second supernatant and an aliquot and quickfreeze both

### **Protein purification**

Thaw supernatants on ice, combine all together and determine volume

Use 200ul of Talon metal affinity resin slurry per 600ul of cell extract

Wash resin 5 times with 3 volumes of binding buffer, spin at lowest speed for 30 seconds for each wash

In cold room, combine resin and cell extract (about 400ul resin to 1.2 ml extract in a 2 mL epi)

Rotate for at least 2 hours to let protein bind to resin

Spin at lowest setting for 30 seconds and collect supernatant S1 (unbound proteins)

Wash the resin of contaminants and unbound/loosely bound proteins

500ul of extraction buffer, rotate 30 minutes

Spin for 1 minute and collect supernatant S2

Repeat wash and spin, collect S3  
Add 250ul of elution buffer with 0.2 M imidazole  
Rotate for 30 minutes, spin and collect S4  
Add 250ul of elution buffer with 0.4M imidazole  
Rotate 30 minutes, spin and save S5  
Add 250ul of elution buffer with 0.6M imidazole  
Rotate 30 minutes, spin and save S6  
Add 250ul of elution buffer with 1.0M imidazole  
Rotate 30 minutes, spin and save S7

### **Xis assay**

Run a protein gel and do a coomassie stain on the four xis containing fractions to determine purity. Combine clean fractions and titrate activity with an excision assay with an array of xis dilutions.

### **Buffers**

#### Binding Buffer

50 mM Tris-HCl pH8  
300 mM KCl  
10% glycerol

#### Elution Buffer

50 mM Tris-HCl pH 7.6  
375 mM KCl  
10% glycerol  
0.2 to 1.0 M Imidazole