



ISSUE 17, 2009

BIOLOGICAL SEPARATIONS

NEWS BULLETIN FOR AXIS-SHIELD DENSITY GRADIENT MEDIA

The risks and problems of using Percoll® for the isolation of cells – comments from publications

Important properties of gradient media

The standard solution of Percoll® that is used for most cell separations is a colloidal silica suspension that is stabilized by polyvinylpyrrolidone (PVP). The solution contains trace amounts of free PVP, which is a known toxin. Percoll® is not endotoxin-tested, but measurements indicate that the levels are at least 6 EU/ml.

In contrast all Axis-Shield density gradient media are produced in facilities under strict FDA and EU cGMP compliance to the EU Pharmacological Standard of <1 EU/ml. Certificates of Analysis are available for every batch and the measured endotoxin levels are routinely <0.13 EU/ml.

The importance of low endotoxin levels is indicated by the following effects of lipopolysaccharide: *In vivo* endotoxin (lipopolysaccharide) causes a variety of inflammatory responses, hypotension, changes to leukocyte populations, intravascular coagulation, shock and death. *In vitro* endotoxin interacts with CD14 and other receptors; causes cytokine production in monocytes and macrophages; activates complement and coagulation cascades and acts as a B cell mitogen.

The wide range of mammalian and non-mammalian cells that demonstrate sensitivity to Percoll® suggests that possible effects on any cell type from this medium may need to be considered. The text below describes some of the specific problems (purity, yield and/or function) that have been encountered with Percoll® gradients.

General effects on mammalian cells

One of the first publications which warned about the possible deleterious effects of Percoll® was published over 25 years ago by Wakefield J.S. et al [1] who observed that peritoneal macrophages, hepatic Kupffer cells and testis Leydig cells phagocytosed the colloidal silica particles. Even though phagocytosis can be reduced at 4°C, the adherence of macrophages to plastic was inhibited by the Percoll® at this temperature and the authors suggested that a number of surface phenomena may be affected by a loose association of the colloidal silica particles to the cell surface.

Lung cells

Viscardi, R.M. et al [2] were concerned that Percoll®, a gradient medium that is commonly used for type II pneumocyte cell isolation, may be toxic to cells. They favoured the use of Nycodenz®, emphasizing that it is a true solution, transparent, not metabolized by cells, and nontoxic to cells.

Gastric parietal cells

Benn et al [3] used continuous gradients favoured the use of Nycodenz® for these cells. Purities >95% were obtained; Percoll® gave much inferior recoveries and purity of parietal cells and also stimulated acid and cAMP secretion by gastric cells. More recently the method has been OptiPrep™, which gave even greater purity [4].

Progenitor cells

Mayanagi et al [5] emphasized the lack of toxicity of Nycodenz®, compared to Percoll®, and the avoidance of positive selection with antibody-beads was important in recovering viable cells.

Percoll is a trademark of GE Healthcare companies

Web:

www.axis-shield-density-gradient-media.com

AXIS-SHIELD

PO Box 6863 Rodelokka

N-0504 Oslo

Norway

Phone: +47 24 05 60 00

Fax: +47 24 05 60 10

Email: bjh@no.axis-shield.com or

john@jgrescon.fsbusiness.co.uk

Erythrocytes

The deformability and rigidity of reticulocytes and erythrocytes may be compromised by the use of Percoll®. There is evidence that the colloidal silica particles adhere to the surface of erythrocytes (even after washing) and cause progressive hemolysis [6].

Bacteria

The recovery, purity and functional integrity of the bacteria isolated from soil in Nycodenz® is superior to isolation in CsCl and Na₂WO₄ and Percoll® [12,13]. Isolation of obligate intracellular bacterium in OptiPrep™ has shown better yields and improved resolution from host cell components than with Percoll® [14].

Sea urchins

For the fractionation of coelomocytes sucrose gradients pose the problem of high osmolality; the silica particles of Percoll® may be engulfed by potentially phagocytic cells and they also tend to adhere to cell surfaces. The cells also have an extreme sensitivity to lipopolysaccharide [15.] Opti-Prep™ posed none of these problems and provided excellent resolution of the different cell types.

Protozoa

Nycodenz® gradients were shown to be superior to those of Percoll® for *Cryptosporidium* oocyst purification, when considering both recovery and particulate load [7]. Iodixanol gradients demonstrate an ability to resolve different cell types that is lacking in other density gradient media in the case of both *Sarcocystis neurona* [8] and *Enterocytozoon bieneusi* [9]. The use of Nycodenz® to enrich for either gametocytes or ookinetes from cultures of *Plasmodium falciparum* greatly improved the viability compared to those purified in Percoll® [10]. Although Percoll® gradients were able to provide a purified sporocyst fraction from *Toxoplasma*, because these particles do not all band in a discrete manner in such gradients, they were unable to provide the simultaneous isolation of a pure oocyst wall fraction possible with OptiPrep™ [11].

Bibliography

1. Wakefield, J.S., Gale, J.S., Berridge, M.V., Jordan, T.W. and Ford, H.C. (1982) *Is Percoll innocuous to cells?* Biochem. J., **202**, 795-797
2. Viscardi, R.M., Ullsperger, S. and Resau, J.H. (1992) *Reproducible isolation of type II pneumocytes from fetal and adult rat lung using Nycodenz density gradients* Exp. Lung Res., **18**, 225-245
3. Benn, S.E., Canfield, S.P. and Curwain, B.P. (1987) *Preparation of purified parietal cells from rat using Nycodenz density gradients* J. Physiol., **391**, 9P
4. Chew, C. S., Parente, J. A., Chen, X., Chaponnier, C. and Cameron R. S. (2000) *The LIM and SH3 domain-containing protein, lasp-1, may link the cAMP signaling pathway with dynamic membrane restructuring activities in ion transporting epithelia* J. Cell Sci., **113**, 2035-2045
5. Mayanagi, T., Kurosawa, R., Ohnuma, K., Ueyama, A., Ito, K. and Takahashi, J. (2003) *Purification of mouse primordial germ cells by Nycodenz* Reproduction, **125**, 667-675
6. Kamiyama M. (2002) *Personal communication*
7. Chesnot, T. and Schwartzbrod, J. (2004) *Quantitative and qualitative comparison of density-based purification methods for detection of Cryptosporidium oocysts in turbid environmental matrices* J. Microbiol. Meth., **58**, 375-386
8. Ellison, S.P., Greiner, E. and Dame, J.B. (2001) *In vitro culture and synchronous release of Sarcocystis neurona merozoites from host cells* Vet. Parasitol., **95**, 251-261
9. Zhang, Q., Singh, I., Sheoran, A., Feng, X., Nunnari, J., Carville, A. and Tzipori, S. (2005) *Production and characterization of monoclonal antibodies against Enterocytozoon bieneusi purified from rhesus macaques* Infect. Immun., **73**, 5166-5172
10. Carter, E.H., Suhrbier, A., Beckers, P.J.A. and Sinden, R.E. (1987) *The in vitro cultivation of P. falciparum ookinetes, and their enrichment on Nycodenz gradients* Parasitology, **95**, 25-30
11. Everson, W. V., Ware, M. W., Dubey, J. P. and Lindquist, H. D. L. (2002) *Isolation of purified oocyst walls and sporocysts from Toxoplasma gondii* J. Eukaryot. Microbiol., **49**, 344-349
12. Rockne, K.J., Liang, W., Young, L.Y. and Taghon, G.L. (2003) *Toxicity of density separation media to Escherichia coli and Myobacterium strain PC01: implications for density-separation of soils and sediments* FEMS Microbiol. Ecol., **43**, 185-189
13. Robe, P., Nalin, R., Capellano, C., Vogel, T.M. and Simonet, P. (2003) *Extraction of DNA from soil* Eur. J. Soil Biol., **39**, 183-190
14. Henriquez, V., Rojas, M.V., and Marshall, S.H. (2003) *An alternative efficient procedure for purification of the obligate intracellular fish bacterial pathogen Piscirickettsia salmonis* Appl. Environ. Microbiol., **69**, 6268-6271
15. Gross, P.S., Clow, L.A., Smith, L.C. (2000) *SpC3, the complement homologue from the purple sea urchin, Strongylocentrotus purpuratus, is expressed in two subpopulations of the phagocytic coelomocytes* Immunogenetics, **51**, 1034-1044