

DENSITY GRADIENT MEDIA

Why should I change to OptiPrep for the purification of virus?

Compared to CsCl and sucrose there are procedural advantages to the use of OptiPrep

- OptiPrep is a sterile solution of 60% (w/v) iodixanol; simply dilute it with saline to prepare sterile gradient solutions. It is the only gradient medium manufactured under strict FDA and EU cGMP compliance.
- CsCl and sucrose are toxic to cells and must be removed prior to infection of cells.
- Iodixanol is non-toxic to cells; it has very low endotoxin levels (<1 EU/ml); measured levels on each batch are usually <0.13 EU/ml. No dialysis is required before re-infection of cells¹.
- CsCl must be removed prior to HPLC or gel electrophoresis; iodixanol rarely needs removing prior to further processing.

CsCl gradients lead to major reductions in infectivity of all viruses

- The recovery of rAAV from iodixanol gradients is at least 10x greater than from CsCl gradients².
- Infectivity titre is up to 10x higher in iodixanol isolates compared to CsCl isolates of rAAV² and adenovirus³.

Sucrose gradients cause changes to the virus surface and poor resolution from cellular components

- Iodixanol gradients give much better resolution of HIV-1⁴ and Herpes virus (HHV-6A)⁵ from membrane vesicles; in iodixanol gradients contamination of HIV-1 by extracellular Vif is also substantially reduced, while infectivity is unimpaired⁴.
- Sucrose gradients lead to loss of surface glycoproteins from retroviruses⁶. Electron microscopy of human endogenous retrovirus particles recovered from iodixanol gradients shows a much better retention of surface architecture⁷. Sucrose gradients strip the majority of glycoprotein gB from Kaposi's sarcoma associated herpes virus; this did not occur in Nycodenz gradients⁸ (iodixanol is a close derivative of Nycodenz).

The Axis-Shield virus database

- Over 35 different viruses have been purified using OptiPrep. Presently 25 Application Sheets giving detailed protocols for virus purification (with full bibliographical back-up) are available from the Axis-Shield website: go to www.axis-shield-density-gradient-media.com and click on "Viruses", then select the virus of interest from the Index.

Web:
www.axis-shield-density-gradient-media.com

AXIS-SHIELD

PO Box 6863 Rodelokka
N-0504 Oslo
Norway

Phone: +47 24 05 60 00

Fax: +47 24 05 60 10

Email: bjh@no.axis-shield.com or
john@jgrescon.fsbusiness.co.uk



Bibliography

1. Krijnse-Locker, J., Kuehn, A., Schleich, S., Rutter, G., Hohenberg, H., Wepf, R. and Griffiths, G. (2000) *Entry of the two infectious forms of vaccinia virus at the plasma membrane is signaling-dependent for the IMV but not the EEV* Mol. Biol. Cell, **11**, 2497-2511
2. Hermens, W.T.J.M.C., Ter Brake, O., Dijkhuizen, P.A., Sonnemans, M.A.F., Grimm, D., Kleinschmidt, J.A. and Verhaagen, J. (1999) *Purification of recombinant adeno-associated virus by iodixanol gradient ultracentrifugation allow rapid and reproducible preparation of vector stocks for gene transfer in the nervous Human Gene Ther.*, **10**, 1885-1891
3. Manninen, A., Verkade, P., Le Jay, S., Torkko, J., Kasper, M., Fullerkrug, J. and Simons, K. (2005) *Caveolin-1 is not essential for biosynthetic apical membrane transport* Mol. Cell Biol., **25**, 10087-10096
4. Dettenhoffer, M. and Yu, X-F. (1999) *Highly purified human immunodeficiency virus type 1 reveals a virtual absence of Vif virions* J. Virol., **73**, 1460-1467
5. Hammarstedt, M., Ahlqvist, J., Jacobson, S., Garoff, H. and Fogdell-Hahn, A. (2007) *Purification of infectious human herpesvirus 6A virions and association of host cell proteins* Virol J. **4**:101
6. Palker, T.J. (1990) *Mapping of epitopes on human T-cell leukemia virus type 1 envelope glycoprotein* In: Human Retroviral HTLV (ed. Blattner, W.A.) Raven Press, NY, pp 435-445
7. Møller-Larsen, A. and Christensen, T. (1998) *Isolation of a retrovirus from multiple sclerosis patients in self-generated iodixanol gradients* J. Virol. Meth., **73**, 151-161
8. Zhu, F.X. and Yuan, Y. (2003) *The ORF45 protein of Kaposi's sarcoma-associated herpesvirus is associated with purified virions* J. Virol., **77**, 4221-4230